



ISOLATION AT THE HIGHEST LEVEL

Leading Isolation Technology

Simply Chemistry. Professional-quality Recrystallization.

Demystify the isolation process with our easy-to-use leading technology. Our isolation systems walk you through the recrystallization process with everything from glassware and vacuum pumps to flow-monitoring control panels.

With our isolation systems, you can produce the purest CBD isolate possible. Whether you're seeking a turnkey manual solution or an automated system that fits seamlessly with the rest of your equipment, Colorado Extraction Systems makes the recrystallization process easier than ever. Our isolation systems are available in a manual 2L version, which comes with our ISO Kit, and an automatic 20L version to perfectly fit your throughput needs.


- Manual and automatic models available
- 2L or 20L capacity
- Processes up to 44 pounds per hour
- User-friendly, intuitive design
- Highly adaptable



ISO2: Manual small batch process utilizing custom designed isolate to recrystallization equipment to maximize yield.

ISO20: Double jacketed Borosilicate glass reactor with Teflon mixer, electronic lift & tilt, 180 degree reactor tilt for ease of operation, heating chilling, and vacuum included.

Model	Capacity (L)	Footprint (ft x ft)	Processing Speed (kg/hr)	Power Requirements
ISO2	2L	4' x 2'	-2 kg	35A 120/1/60
ISO20	20L	8' x 3'	3-20 kg	60A 220/1/60



ISO2



ISO20

Pure Processing Power™

SCHEDULE A DEMO

coloradoextraction.com
800.513.5560 | sales@coloradoextraction.com

